**User instruction.**

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**EN 1.** Ensure that the Velcro Attachment **A1** on the arm part is closed according to markup given by the fitter and that the chest strap **B** and shoulder strap **C** are easily attached to each other. In order to reach the shoulder strap when stretched. Be placed in front of a mirror when putting on the orthosis.

**2.** Pull on the arm part **A** from below and pull it over the shoulder. One should feel that the upper arm follows together a bit up.

**3.** Grip the chest strap **B** and make sure it is placed right over your back. Attach the chest strap to the shoulder part **D** with the black marks against each other. At this stage, the breast band **B** should be stretched in its full length.

**4.** Tighten/adjust the chest strap so that it feels comfortable with the help from the Velcro fastener at the plastic buckle **B2**.

**5.** Grip the shoulder strap **C**. Make sure the shoulder strap is loosened from the chest strap. Slowly pull the shoulder strap outwards/downwards until the desired support and "lift" is reached in the arm/shoulder. Then forward it so that the middle Velcro patch **C1** attaches to the armpit. Finally attach the shoulder strap's outer Velcro fastener **C2** at the front. If necessary, the pillow **E** can be attached under the arm.

**6.** Done. For further information or other questions, contact your fitter.

**NordiCare Ortopedi & Rehab AB**  
Solrosvägen 1, SE-263 62 Viken, Sweden  
Phone +46 42 35 22 20, Fax +46 42 35 22 21  
info@nordicare.se, www.nordicare.se

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**EN** If a stinging sensation, rash or other changes in skin color occur, discontinue use and contact your doctor or other medically trained person. Only medically trained persons are to prescribe and recommend the use of the product. For best results, follow instructions for use. Periodically check the condition of the Velcro to ensure they are free from grease and damage. Use in combination with other medical products should be supervised by qualified medical practitioner. Serious incidents shall be reported to NordiCare and the competent authority.

For more information please consult our product catalogue or visit [www.nordicare.se](http://www.nordicare.se). For further questions, contact your doctor, health practitioner or medical equivalent - or NordiCare: +46 42-35 22 20, info@nordicare.se

This packaging can be recycled as paper and the bag as plastic.

DISPOSAL OF THE PRODUCT  
Do not dispose the product in the nature. Dispose of the product in accordance with the waste laws of your municipality and/or country. The product does not present a risk of contamination.

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**EN** HemiSafe is a functional shoulder orthosis specifically designed for use when in sublaxation or shoulder pain after stroke.

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**7.** Valmis.

#### EN

**1.** To adjust the arm part **A** to the correct size and function, close the arm part below the elbow and pull it gently upwards. Change it to reach the desired fit. The arm section should be tight but not too tight. Please mark the placement **A1** with a pen to help the user. Pull the arm section from the bottom to ensure that the upper arm is "lifted" up in the appropriate position.

**2.** Adjust the chest band **B** which function is to hold the orthosis in place. Attach the end of the chest strap loosely to the shoulder part **D** using the black markings against each other.

**3.** Before tightening the chest strap, check the length: open the D-ring bracket **B1** with crocodile gap and shorten the chest strap if it is too long. The position of the D-ring should be just in front of the armpit.

**4.** Then tighten the front band of the chest at **B2** so it is comfortable around the chest, not too tight.

**5.** Tighten the shoulder strap **C** smoothly and gently outwards/downwards so that the desired support and "lift" is achieved in the arm/shoulder. The shoulder strap function is to fix the arm/shoulder to the correct position. Make sure the middle Velcro attachment **C1** on the shoulder strap attaches to the chest strap. Then attach the shoulder strap outer part **C2** to the front.

**6.** If necessary, the pillow **E** can be attached under the arm when increased abduction is desired.

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**EN** **1.** Circumference in cm around the upper arm measurement is to be taken where the upper arm is the largest.  
**2.** The circumference measure around the chest, chest size measurements just below the armpit.

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**EN** Shoulder joint orthosis after stroke to limit sublutation of the shoulder joint.

**Contraindications:** Use caution in case of reduced sensation.

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### Material Composition

Polyester 62%, Polyamide 36%, Elastan 2%. Latex free.