



QuickSling Arm sling Art no. 20118-00-50

User Instructions

Double arm sling

1. Cut to desired length.
2. Attach Velcro to both ends of the strap.
3. Make a loop for application near the elbow.
4. Pass the sling behind the back over the opposite shoulder.
5. Wrap the sling around the wrist and adjust the position of the hand by attaching the Velcro to the sling.

Indications: For relieving shoulders/arms in the event of acute injury/trauma.

Contraindications: For patients with sinus carotid sensitivity avoid pressure from the orthosis in the carotid sinus area.

Material:

53% Polyurethane, 38% Polyester, 9% Nylon.
Latex free.

