Ballerina

art no 50103

POSTURE VEST


DISPOSAL OF THE PRODUCT
Do not dispose the product in the nature. Dispose of the product in accordance with the waste laws of your municipality and/or country. The product does not present a risk of contamination.

This packaging can be recycled as paper and the bag as plastic.

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
EN Functional posture vest that improves
posture and ease of respiration. Used to
assist with cases of osteoporosis, increased
thoracic kyphosis, compressed vertebrae,
and clavicle fractures.

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EN Make sure that the shoulder straps are crossed through the loop **A** in the back and through the loops **B** on the waistband.

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EN 1. Put on the posture bandage over the shoulders like a waistcoat.

2. Tighten the waistband and secure the Velcro (figure 1). To avoid pressure on the diaphragm, fasten them obliquely downwards.

3. Take hold of the adjustable straps and pull them out and forwards until the shoulders are pulled back to the appropriate position, (figure 2). Secure the Velcro fastening at the front, (figure 3).

4. If the adjustable straps are too long, open the y-gap, and cut the strap to the required length, (figure 5, C).

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Washing Instructions

Machine wash, 40°C (104°F). Close the Velcro. Use a washing net. Dry flat. Do not bleach, not dry clean. Do not iron or tumble dry.

Material

85% Polyester, 7% Spandex, 6% Nylon, 2% Cotton. Latex free.

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EN Functional thoracic- and shoulder posture bandage. Thoracic kyphosis and posture related problems.

Contraindications: Unstable thoracic fracture.

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EN If a stinging sensation, rash or other changes in skin color occur, discontinue use and contact your doctor or other medically trained person. Only medically trained persons are to prescribe and recommend the use of the product. For best results, follow instructions for use. Periodically check the condition of the Velcro to ensure they are free from grease and damage. Serious incidents shall be reported to NordiCare and the competent authority.

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For more information please consult our product catalogue or visit www.nordicare.se. For further questions, contact your doctor, health practitioner or medical equivalent - or NordiCare: +46 42-35 22 20, info@nordicare.se.