











## KRB - Telescopic Knee ROM

EN	<b>KNEE ORTHOSIS</b>
DE	<b>KNIEORTHESE</b>
FR	<b>ORTHÈSE DU GENOU</b>
ES	<b>ORTESIS DE RODILLA</b>
IT	<b>ORTESI PER GINOCCHIO</b>
NL	<b>KNIEORTHESE</b>
DK	<b>KNÆORTOSE</b>
FI	<b>POLVIORTOOSI</b>
SE	<b>KNÅORTOS</b>
NO	<b>ORTOSE, KNE, ORDINÆR</b>
PL	<b>ORTEZA STAWU KOLANOWEGO</b>
GR	<b>ΝΑΡΘΗΚΑΣ ΓΟΝΑΤΟΣ</b>

<b>1</b>	
<b>2</b>	
<b>3</b>	
<b>4</b>	
<b>5</b>	
<b>6</b>	
<b>7</b>	
<b>8</b>	
<b>9</b>	
<b>10</b>	
<b>11</b>	
<b>12</b>	
<b>13</b>	
<b>14</b>	

#### ENGLISH

##### Product Description

A post op knee brace to provide controlled range of motion at the knee joint. Telescoping bars can be adjusted to customise brace length.

##### Usage Information (See diagram)

- Unclip plastic buckles and lay brace flat. Note the hook strips at positions A & B. Unhook to allow adjustments of straps, ensuring final buckle position is over the pad for comfort.
- Position the brace with the knee centered between hinges. Press buttons C & D to extend telescoping hinge bars to accommodate leg length. Ensure both hinges are set to the same length above and below knee.
- Position the hinge bars medially and laterally to the leg and centre hinge mechanism at the knee point as shown.
- Fasten the two strap/buckles closest to the knee.
- Loosely fasten the strap back around the brace, using the hook/loop "Y" tabs. (Length can be cut and adjusted later if required)
- Fasten the remaining two buckles and loosely fasten the straps, in the same way as step 5.
- Now tighten all the straps, maintaining the medial and lateral positions of the hinge bars.
- Removed hook/loop "Y" tabs and trim straps if required. Take care to leave enough length for future adjustments if required.
- Extension limit settings may be selected between 0 and 90 degrees by pressing the button and sliding to desired position.
- Flexion limit settings may be selected between 0 and 120 degrees.
- The hinge may be locked by sliding the drop lock button into the locked position. At any one of four positions 0, 15, 30, 45 degrees of flexion.
- Plastic ties can be used to avoid the flexion/extension settings being adjusted. Slide one of the plastic ties provided under side of the button shown.
- Bend plastic tie around and underneath opposite side of the button, ensuring it is secure.
- Tighten and cut off excess.

##### Indications for Use

- Following surgery of the knee joint.
- For immobilisation after knee joint injuries and subsequent mobilisation.
- Knee ligament injuries.

##### Contraindications

- Do not use over open wounds.
- Any known allergy and/or hypersensitivity to any material listed in product composition.

##### Warnings and Precautions

- We recommend the initial fitting of this brace be conducted by a suitably qualified healthcare professional who will advise the period of use.
- Carefully read all instructions and warnings prior to use.
  - Follow all instructions to ensure proper performance of the brace.
  - Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area.
  - Do not re-use for another patient, doing so risks cross infection and can compromise product integrity.
  - Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
  - Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
  - The durability of the brace may be compromised by certain factors, e.g., objects with sharp edges or damage to the hook and loop fasteners. To prevent this the hook and loop fasteners should always be fastened when the brace is not being worn or when being washed.

##### Washing and Care Instructions

- Close over all Hook/Loop straps.
- Hand wash pads/straps in cool water using mild detergent.
- Rinse and dry flat.
- Do not spin or tumble dry.
- Do not use fabric conditioner at any time.
- Check for any degradation that may affect use.
- Periodically check Hook/Loop fasteners are operating securely.

##### Product Family Composition

Aluminium, Nylon, Fibreglass, Low carbon steel, PU, Polyester, POM.

##### Storage and Transport Conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

##### Recycling and Disposal

Packaging and constituent parts should be recycled or safely disposed in accordance with local or national laws.

##### Serious Incident

Report any serious incident to the manufacturer and the competent authority of the EU Member state, or country in which you reside.

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©2023 Promedics Orthopaedics Ltd. All Rights Reserved.