



## Carpal Orthosis Grade 3

Stabilizes the carpal joint with extra stable support. Post-op and conservative treatment for medium to severe sprains and medium arthrosis. Provides rigid lateral and medial support and prevents carpal joint movement.

The orthosis uses a thermoplastic stabilizer system. By removing the thermoplastic sheet, you reduce from grade 3 to grade I. The orthosis comes with a waterproof and breathable cover made in ScandiTex.

### Indications

- Protects against re-injury
- Joint position sense proprioception
- Reduces the carpal joint deviations and hyperextension
- Post-operatively
- Medium to severe sprains
- Medium arthrosis

### Size

Stocked in five Sizes: XS, S, M, L and XL  
Available for both left and right carpal  
Contact Scandi Orthopedic if other size is needed

### Material

Neoprene, aquaplastic, velcro

### Washing instructions

Remove the thermoplastic sheet, as it cannot be washed in heated water. Clean the orthosis with mild detergent, by hand, maximum 30 °C.

# Size calculation


## Carpal Orthosis Grade 3



### Measurement

Take measurements in centimeters with the animal standing.

1. Measure the circumference of the carpal joint
2. Measure the distance from the elbow to the top of the paw



### Size guide

Select the appropriate size below.

Size	1. (cm)	2. (cm)
XS	Up to 10	6.5 – 11
S	Up to 12	9 – 16
M	Up to 15,5	10.5 – 20
L	Up to 19,5	13 – 25
XL	Up to 22,5	15 – 30