





Medical Protection Cover Front Leg with ScandiTex

Professional alternative to a collar. Protects the dog's front leg so that the dog cannot bite or lick on a surgical wound, drainage, stitches and itches. Made with ScandiTex, a waterproof, strechy and breathable fabric. No irritating seams on the inside.

Indications

- Instead of stressful collar
- Post-op.
- · Skin diseases, irritation and wounds
- Drainage
- Over bandages
- Can be used 24/7

Size

Stocked in seven sizes: XXXS, XXS, XS, S, M, L and XL. Contact Scandi Orthopedic if other size is needed.

Material

ScandiTex, Velcro

Washing instructions

Clean the product with mild detergent, by hand, maximum 40 $^{\circ}$ C.

Scandi Orthopedic Ringögatan 5 417 07 Göteborg, Sverige Tel: (+46) 031-214 128 E-mail: info@scandiop.se www.scandiop.se

Size calculation

Medical Protection Cover Front Leg with ScandiTex



Measurement

Perform the measurement with the dog standing. Use the metric system.

- I. Measure the circumference of the chest at its widest (just behind the front legs).
- 2. Measure the distance between the armpit and the top of the paw.



Sizing

Choose an appropriate size below.

Size	I. (cm)	2. (cm)
XXXS	25,5 - 33,5	10,5 - 13
XXS	38,5 - 53,5	13,5 - 19,5
XS	37,5 - 45,5	20 - 24,5
S	52 - 67	23,5 - 29,5
M	65 - 88	27 - 33,5
L	76,5 - 107,5	31 - 37
XL	81 - 117	34 - 42

Scandi Orthopedic Ringögatan 5 417 07 Göteborg, Sverige Tel: (+46) 031-214 128 E-mail: info@scandiop.se www.scandiop.se