

Your silicone for any application

Indication	Dynax® clear	Dynax® heavy body	Dynax® light	Dynax® mono	Dynax® putty	Zerosil® soft
Crown and bridge impressions						
Double mixing technique	•	•	•	•	•	
Correction	•	•	•	•	•	
Sandwich	•	•	•	•	•	
Monophase	•	•	•	•	•	
Inlay and onlay impressions						
Double mixing technique	•	•	•	•	•	
Correction	•	•	•	•	•	
Sandwich	•	•	•	•	•	
Monophase	•	•	•	•	•	
Pick-up impressions e. g. for implants						
Correction	•	•		•		
Monophase	•	•		•		
Functional impressions						
Function				•		
Fixation impressions						
Monophase		•		•		
Relining impressions						
Monophase			•	•		
Pressure spot indicator						
Monophase			•	•		
Impressions for temporaries						
Monophase				•		•
Situation				•		•
Situation impressions						
Monophase						•
Situation						•
Impressions of the opposite jaw						
Monophase						•



Silicones need more than 30 minutes to reach their perfect characteristics. Please make sure that you do not fill the silicone impression immediately after taking it out of the patient's mouth.



Situation impression with Zerosil® soft