

StomiFix - Stoma binder

ART.NO 50217

Soft fixation binder in seamless knitted material. For fixation of ostomy pouches. StomiFix has a tubular shape and consists of two layers. The inner layer has a pouch hole and the outer layer provide support and has a pocket to fixate and hold the bag safely in place, while reducing any sounds. The stoma binder is put on just like regular underwear. It is easy to use and is not visible under clothing. StomiFix can be used every day, night and day, also during sporting activities. Wear your regular underwear underneath or on top of the StomiFix. Suitable for stoma on either right or left side. StomiFix is available in 3-pack in black or white. Height in front 20 cm, in back 20 cm.



Levels of compression: Soft, 1 out of 3.

Indications: To be used for stoma to fixate/keep the ostomy pouch in place.

Contraindications: Ostomy hernia.

Material: 90% Polyamide, 10% Elastane. Latex free.

Washing instructions: Machine wash in 60°C (140°F). The product lasts longer if washed in 40°C (104°F). Possible to tumble dry at low heat. Do not iron. Do not use bleach, chlorine bleach or fabric softener as this could impair the function of the material. Do not dry clean.



Measurements: In order to determine the correct size, measure the hip circumference at its widest point.

| Size | Hip circumference cm | Art.no StomiFix, white, 3-pcs. | Art.no StomiFix, black, 3-pcs. |
|------|----------------------|--------------------------------|--------------------------------|
| XS/S | 65 - 90 | 502178930 | 502178030 |
| S/M | 85 - 110 | 502178940 | 502178040 |
| M/L | 105 - 125 | 502178950 | 502178050 |
| L/XL | 120 - 140 | 502178960 | 502178060 |

