

# Frejka - Hipabduction orthosis


ART.NO 60232

The Hipabduction rail is used to position the hip joints of children having instability in the hip joint. Frejka consists of a flexible polyethylene rail with padded edges. Outside the rail pants of textile.

**Indications:** Congenital hip dislocation/subluxation. Hip dysplasia. **Contraindications:** Dislocated hip.

**Material Composition:** Polyethylene 55%, Polyester 24%, Cotton 17%, Polyamide 4%. Latex free.

**Washing Instructions:** Close all Velcro straps before washing. The Pant is machine washable at 40°C. Wipe the plastic part with a soft damp cloth, use mild detergent.



Size	Width, plastic part (cm)	Measurements between the knees (cm)
X Small	14	15-18
Small	16	17-20
Medium	18	19-22
Large	20	21-24
X Large	22	23-26
2X Large	24	25-28



MD Class I CE