




# Ventilator Circuit

Thank you for purchasing this Armstrong Medical Ventilator Circuit  
Before clinical use, please observe the following guidelines:

## Directions for use:

1. Remove product from packaging ensuring that no separate components become mislaid
2. Before use, verify flow patency and test for leakage
3. This circuit assembly is warranted for use on 1 (one) patient for a maximum of 7-days from date of initial use or removal from primary packaging
4. Contact Armstrong Medical Ltd, or our representative in your country if you are not satisfied with the quality of this product

Key to Symbols:

	Patient End
	Detachable Connection
	adjust length if required

**Keep Package Sealed Until Use**  
**Risk of Blockage by Small Objects**

This product has been manufactured in compliance with relevant parts of ISO 5367 and BS EN ISO 5356-1



Wattstown Business Park  
Newbridge Road  
Coleraine BT52 1BS  
Northern Ireland

Tel: +44 (0)28 7035 6029  
Fax: +44 (0)28 7035 6875  
e-mail: info@armstrongmedical.net



Contains silver phosphate glass antimicrobial technology to protect against microbial growth




Manufactured in the European Union


## Ventilator Circuit

REF AMVC1792/141

For Non-Invasive Ventilation



**Mask Size**  
**Adult**



Drawing not to scale

LOT  


300605  
2010-05