




Ventilator Circuit

Thank you for purchasing this Armstrong Medical Ventilator Circuit
Before clinical use, please observe the following guidelines:

Directions for use:

1. Remove product from packaging ensuring that no separate components become mislaid
2. Before use, verify flow patency and test for leakage
3. This circuit assembly is warranted for use on 1 (one) patient for a maximum of 7-days from date of initial use or removal from primary packaging
4. Contact Armstrong Medical Ltd, or our representative in your country if you are not satisfied with the quality of this product

Key to Symbols:

	Patient End
	Detachable Connection
	adjust length if required

Keep Package Sealed Until Use
Risk of Blockage by Small Objects

This product has been manufactured in compliance with relevant parts of ISO 5367 and BS EN ISO 5356-1



Wattstown Business Park
Newbridge Road
Coleraine BT52 1BS
Northern Ireland

Tel: +44 (0)28 7035 6029
Fax: +44 (0)28 7035 6875
e-mail: info@armstrongmedical.net



Contains silver phosphate glass antimicrobial technology to protect against microbial growth




Manufactured in the European Union


Ventilator Circuit

REF **AMVC1792/145**

For Non-Invasive Ventilation



Mask Size
Small



Drawing not to scale

LOT


300605
2010-05