

INSTRUCTIONS FOR USE AND MAINTENANCE

ELASTIC NECK LOOP FOR SRT

REF.

FTC/TSRT

FEBRUARY 2017 VERSION (TRANSLATED)

1 PRODUCT DESCRIPTION

This elastic neck loop will hold the chest ascender upright and taut during SRT ascent (single rope). The length can be adapted to any stature. It is easy to use, no need to take off the helmet. Practical, compact and ultra-light, it is also quickly stored away.

2 TECHNICAL DATA

The length is adapted with the plastic slider. Flat rubber band 28 mm large, woven from polyester and cotton. The black side is breathable, thus comfortable to the climber.



FTC
PLAY WITH GRAVITY

1355, chemin de Malombre • Z.I. Les Plaines
26780 MALATAVERNE • FRANCE • (+33) 475 528 640
contact@ftc-tree.com • www.ftc-tree.com

3 SETTING UP



Slide the carabiner in the smaller loop, then lay the rubber band flat in front of you (still open), printed side up.



Pass the band behind your neck, laying both carabiner and larger loop on your chest.



Slide the carabiner in the larger loop, overlapping on the smaller one.



Clip the carabiner on the chest ascender.



Adapt the length with the slider.