

Miriam, postoperative bra for fixation



NordiCare[®] +

Miriam, postoperative bra for fixation

Seamless and adjustable bra for securing dressings


- Designed for use after breast surgery/reconstruction and radiotherapy.
- Adjustable and effective fixation of dressings.
- Fully openable at the front and at the shoulder straps for easy application directly after surgery and during the healing process.
- Extra wide across the back and high under the arms.
- Pressure-free cups.



Wide, adjustable shoulder straps with front adjustment. Available in white and black.



Wide racerback and high sides. Soft, seamless under the bust to minimise irritation.



Fully openable with hook-and-eye fastenings at the shoulder straps and center front.



Extra-stretch, pressure-free cups.

Miriam post-op bra

cm	60	65	70	75	80	85	90	95	100	105	110	115
Cup size	A											
	B											
	C			S	M	L	XL	2XL				
	D	XS										
	E											
	F											
	G											

Model	Size	Item No.
Miriam white	XS-2XL	31154-09
Miriam black	XS-2XL	31154-00



Latex-free

MD Class I



Helena Peters
+46 (0)42-450 88 51
helena@nordicare.se