


Technical specification sheet

Product name	: Sensor Sleeve #2
Sheet #	: TECHS02SLEEVE
Colour	: Semi transparant
Artical code	: S02
General purpose	: Hygiene sleeve for sensor
Material	: Vinyl, textured
Technics	: Sealing,
Certificates	: CE : NL/CA01/2006-02593, FDA: Product code: EGZ Classification name: Section 872 1905 Dental X-Ray Film holder Listing number: E591672
Recycle	: No, disposable, one time use only
Supplied as	: Part of introduction kits and refill's



Application

Preventing transfer contamination by body fluids and bacteria from one patient to another, by using sensor and sensor holder. Without the need of disinfecting the sensor itself.

Instructions for using the sleeve

Slide sensor into sleeve, until the end.
Place the sensor/sleeve into the sensorholder

General hygiens notes

The sensor sleeve must be used on one patient and after use be handeld as non recycle waste.

General terms

DenBite® X-ray sensor sleeve are manufactured by DuPhaMed B.V. Patent pending.
© 2005 DuPhamed B.V. - Alphen aan den Rijn, The Netherlands.
All rights reserved, Other Brands and products are Trademarks of their respective holders
All technical specification's are subject to be changed without prior notification.